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Chest and pelvis seatbelt set – Mounting instructions

Description

HP VELOTECHNIK chest and pelvis seatbelts help to prevent falling-out-of-the-seat of persons with limited mobility. The seatbelts can be mounted to ErgoMesh meshseats and Gekko meshseats.

Content

- chest seatbelt or pelvis seatbelt with mounting clamps, preassembled

Assembly

- 1) Locate the position of the seatbelt in height on the back of the seat frame.
- 2) Open the seatmesh strap(s) on the seat's backside in the area of the seatbelt mounting position.

- 3) Dismount the clamp from the seatbelt by removing the loosely assembled screw.

- 4) Slide the clamp around the seat's side tube. The angled tongue of the clamp is the forward one and points inwards.

- 5) Pull together both tongues to have both holes nearly one upon the other. Use pliers if necessary.

- 6) Reattach the seatbelt to the clamp. Don't tighten the screw yet. Check the seatbelt's position in height and alter the position of the cushion part on the perforated plastic belt if necessary.

To do so, loosen the Phillips head screw and pull the belt from the plastic noses. Slide the belt to the desired position, press it back on the noses and screw the Phillips head screw back in.

7) To improve the orientation of the seatbelt the plastic belt shall be pre-bent. First, locate the area of the bending. Unscrew the belt from the clamp. Using a hot air gun heat the area of bending slowly and evenly from both sides. As soon as the plastic reached its deformation temperature bend the belt more than 90°. Hold this position and let the belt cool down.

Caution! The belt is very hot after heating. Use protective gloves to avoid burnings.

8) Reattach the belt to the clamp. Tighten the screw with tightening a torque of 5 – 7 Nm.

9) To mount the seatbelt on the right side of the seat, proceed analogically.

In case there is an accumulator mounted to the seat frame that interferes with the desired seatbelt position, use one of the accumulator adapter mounting screws to fix the seatbelt.

10) Reattach the seatmesh strap to its lock and tighten the strap.

The chest or pelvis seatbelt is ready for use now.