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Mounting instructions for Walking Aid Mount on BodyLink®-seats

Content

- (1) Holder with locking sheet (1b)
- (2) Mounting sheet for holder on the seat
- (3) Quiver
- (4) Mounting plate
- (5) Mounting base
- (6) Quiver adapter sheet
- (8) Bolts
 - o Socket head bolt M4x18mm, ISO7380
 - o Round head bolt M5x14, ISO7380
 - o Round head bolt M5x20, ISO7380
 - o Round head bolt M5x25, ISO7380
 - o Round head bolt M5x35, ISO7380
 - o Socket head bolt M8x20, DIN7984
 - o Socket head bolt M8x40, DIN7984
 - o Socket head bolt M8x45, DIN7984
- (9) Self locking nut DIN 985 M5
- (10) Washers
 - o Large washer for M4, DIN9021
 - o Large washer for M5, DIN9021
 - o Washer for M5, DIN125
- (11) Adhesive sticker

This Walking Aid Mount can be used for carrying walking aids on Scorpion-tricycles by HP VELOTECHNIK.

It is possible to use one Walking Aid Mount on each side of a tricycle with BodyLink®- seats.

For ErgoMesh®- and Gekko-seats there is one special kit of walking aid mounts with the part no. GHEM available.

It is possible to handle the system with only one hand: The stand of the walking aid is stuck in the quiver while the holder fixes the shaft with diameters from 20 to 30 mm. The spring steel locking sheet prevents the holder from accidentally opening and loosing the walking aid while riding the tricycle.

The mounting situation of the Walking Aid Mount varies depending on the configuration of the tricycle -depending on depending on the type and individual configuration.

Hint: Before fixing all screws tightly apply both, quiver and holder to the bike. After adjusting the parts to your needs, tighten all screws properly.


I: Mounting of the quiver adapter sheet

If your tricycle is equipped with an electric assist system, it may be necessary to mount the quiver (3) with the adapter sheet (6) on the battery adapter sheet (14).

- Take off the battery from the battery rail
- Disassemble the battery rail from the battery adapter sheet (14).
- Use the large washer for M5 (10) with the M5x16 socket head bolt (8). Put it through the

rear end of the slot in the battery adapter sheet (14).

- **Alternatively** you can attach the quiver adaptor sheet (6) directly to one of the M6 bolts that holds the battery rail. Then it will be necessary to re-drill the hole in the long end of the sheet (6) with a 6 mm drill.
- Attach the quiver adapter sheet (6) with the long end to the battery adapter sheet (14). A M5 washer (10) and a self-locking nut (9) are used to fix it on the bolt.


2: Mounting of the quiver

A: Attaching the quiver (3) on the quiver adapter sheet (6)

The procedure of mounting the quiver on the adapter sheet is basically the same on Gekko and Scorpion-tricycles. See picture below.


- Use the M5x14 round head bolt (8) with a large washer for M5 (10). Put both through the mounting plate (4) into the counterbore.
- Put the mounting plate (4) with the bolt (8) and the washer (10) into the counterbore of the quiver (3).

- Put the bolt (8) through the hole in the quiver adaptor sheet (6), apply one washer for M5 (10) and fix it with a self locking nut M5 (9). Take care that everything is centred in the hole of the quiver (6)

B: Scorpion tricycles: Mounting the quiver on the rear swing arm

Follow this procedure when you have no rear rack or battery mounted.


- On Scorpion, Scorpion fs 20 and Scorpion fx 20 remove the M8 round head bolt from the swing arm axle (13) on the side where the quiver (3) will be mounted.
- Use the socket head bolt M8x20 and put it into the mounting plate (4).
- Put the mounting plate (4) with the bolt (8) and into the counterbore of the quiver (3).
- On the other side, place the mounting base (5) with the shouldered flange into the hole of the quiver.
- Put the bolt (8) into the axle of the rear swing arm on Scorpion, Scorpion fx and Scorpion fs 20 – respectively the swing arm bolt (12) on Scorpion fs 26 or Scorpion plus tricycles. Take care that the mounting plate (4) and the mounting base (5) are centred in the hole of the quiver (3).
- In the end when all parts are mounted and aligned, **remember to tighten the bolt (8) with a torque of 17 – 19 Nm!**

Mounting instructions for Walking Aid Mount on BodyLink®-seats

C: Scorpion tricycles: Mounting the quiver with the rear rack

Follow this procedure when you have a rear rack and no battery mounted.

Attaching the quiver on the swing arm axle with a rear rack mounted is basically the same like described in "B" with no rear rack.

- First you have to release the countersunk bolt from the rear rack (11) and the swing arm bolt (12) / swing arm axle.
- Put the bolt M8x45 (8) through the mounting plate (4), the quiver (3) and the mounting base (5). Then put the screw with all the parts through the hole of the rear rack (11) into the swing arm bolt or (depending on the tricycle model you ride) the swing arm axle.
- In the end when all parts are mounted and aligned, **remember to tighten the bolt (8) with a torque of 17 – 19 Nm.**

3: Mounting of the holder

For mounting the holder (1) on the BodyLink®-seat use the mounting sheet (2). This sheet can be used for attaching the holder on the right or left side of the tricycle.

- Possibly you will have to drill through the holes M6 for the headrest mount with a 5 mm drill to open the holes on the back side of the seat.

- Depending on the accessories you use, like headrest, rear light mount or bottle holder, you will need different length of the ISO 7380 round head bolts (8) M5x20, M5x25, M5x35.
- Put the large washer (10) on the bolts and then stick them through the holes in the back rest of the seat.
- Put the adapter sheet for the holder (2) on the screws and fix it with the self locking nuts (9) for M5 using M5 washers (10).
- If you have a rear light mounted on the seat you can fix the adapter sheet (2) by re-attaching the rear light mount.
- Then we recommend to use thread locker to secure the connection.

- For mounting of the holder (1), use the M5 round head bolt (8) with a large washer for M5 (10). Put them into the counterbore of the mounting plate (4).

- Put the mounting plate (4) into the holder (1) and insert the bolt into the hole of the mounting sheet (2). The holder faces to the rear.
- Attach a washer (10) and a self locking nut (9) for M5 on the bolt to fix the holder (1) after aligning it horizontally.

4: Handling of the Walking Aid Mount

First you have to open the locking sheet sheet (1b) on the holder (1). Put the stand of the walking aid into the quiver (3) and then push the shaft into the rubber part of the holder. The locking sheet (1b) swivels around by itself and the walking aid will remain in its position.

Now push the locking sheet (1b) until it clicks into place. Only if the locking sheet (1b) is locked in the secure position, the walking aid is protected from being released accidentally!

Caution: Only ride your bicycle when the securing sheet is locked in secure position to prevent the walking aids from being released accidentally out of the holder.

To release the walking aid, pull the locking sheet slightly on the black latch. Do not open the holder completely. By taking out the walking aid, the securing sheet swivels around and the holder opens completely.

Attention: Depending on the mounting situation and the mounting side, it might be necessary to turn the quiver away when folding the tricycle to avoid collision with the frame or other parts.

Caution: Make sure the safety of your bicycle will not be affected by the Walking Aid Mount and by taking walking aids with you.